

# Ensuring Quality Data for Organochlorine Pesticide Analysis by Using an Optimized SPE Cartridge for Sample Extract Cleanup

Jason Hoisington, Jason Thomas, Alexis Shelow, Colton Myers

Restek Corporation

NEMC 2025

# EPA 8081 Pesticides Analysis Troubles

## Problems resulting from complex sample matrices

- Quantitation inaccuracies
  - Calibration interval reduction
  - Frequent inlet degradation check failures
- 
- Increased maintenance requirements
  - Increased downtime
  - Reduced throughput
  - Diminishing profitability

# EPA 8081 Pesticides Analysis Problems

## Problems resulting from complex sample matrices

### ➤ Chromatographic interferences

- Coelution from coextracted matrix compounds – e.g., chlorophenols
- Massive chromatographic interferences – e.g., sulfur
- Carryover/ghost peaks/crossover contamination – high boilers

### ➤ System inertness effects

- Continuing calibration standards failure
- Inlet degradation check failure

# EPA 8081 Pesticides Analysis Problems

## Problems resulting from complex sample matrices

### ➤ Chromatographic Interferences


- **Coelution from coextracted matrix compounds – e.g., chlorophenols**  
Methods 3620/3610/3630
- **Massive chromatographic interferences – e.g., sulfur**  
Method 3660
- **Carryover/ghost peaks/crossover contamination - high boilers**  
Method 3640

# EPA 8081 Pesticides Analysis Problems

## Problems resulting from complex sample matrices

### ➤ Chromatographic Interferences

- Coelution from coextracted matrix compounds – e.g., chlorophenols




# EPA 8081 Pesticides Analysis Problems

## Problems resulting from complex sample matrices

### ➤ Chromatographic Interferences

- Massive chromatographic interferences – e.g., sulfur

Method 3660




# EPA 8081 Pesticides Analysis Problems

## Problems resulting from complex sample matrices

### ➤ Chromatographic Interferences

- Massive chromatographic interferences – e.g., sulfur



# EPA 8081 Pesticides Analysis Problems




## Problems resulting from complex sample matrices

### ➤ Chromatographic Interferences

- Carryover/ghost peaks/crossover contamination - High Boilers

Method 3640

GPC costly and cumbersome process



Method 3640A: Gel-Permeation Cleanup, part of Test Methods for Evaluating Solid Waste, Physical/Chemical Methods



# EPA 8081 Pesticides Analysis Problems

## Problems resulting from complex sample matrices

- **System inertness effects**
  - **Continuing calibration standards failure**
  - **Inlet degradation check failure**

9.3.3 DDT and endrin are easily degraded in the injection port. Breakdown occurs when the injection port liner is contaminated with high boiling residue from sample injection or when the injector contains metal fittings. Check for degradation problems by injecting a standard containing only 4,4'-DDT and endrin. Presence of 4,4'-DDE, 4,4'-DDD, endrin ketone or endrin indicates breakdown. If degradation of either DDT or endrin exceeds 15%, take corrective action before proceeding with calibration. Unless otherwise specified in an approved project plan, this test should be performed even when DDT and endrin are not target analytes for a given project, as a test of the inertness of the analytical system.

# EPA 8081 Pesticides Analysis Problems

## Problems resulting from complex sample matrices

- **Inlet degradation check failure**
  - **Endrin – metal activity, lack of deactivation**
  - **DDT – accretion of high molecular weight material in inlet and head of column, ex: humic/fulvic acids**

# EPA 8081 Pesticides Extract Cleanup

- **GPC costly and cumbersome process**
  - **Florisil only, no HMW**
    - chlorophenols
  - **CPP only, no TCP**



# EPA 8081 Pesticides Extract Cleanup

## ➤ Combined orthogonal sorbents

- no cumulative effects of combined strong retention

## ➤ GCB with tempered affinity for planars

- High variation and high retention (HCB) reigned in


## ➤ Florisil with moisture control

- Commonly encountered high moisture levels causing > 5% TCP @ 10 mls elution

# EPA 8081 Pesticides Extract Cleanup


Visual difference in cleanup efficacy

(note: colorless does not always mean it's free from activity problems, nor does a colored extract guarantee problems)



# EPA 8081 Pesticides Extract Cleanup

Comparison of DDT breakdown check standard results checked successively after repeated sets of 10 injections of a soil extract cleaned by Florisil only cartridge versus that of GCB stacked with Florisil






# EPA 8081 Pesticides Extract Cleanup

Notice easy  
pass for  
trichlorophenol  
breakthrough  
requirement

| Compound           | Concentration (ng/mL) | Average Recovery (%) | %RSD |
|--------------------|-----------------------|----------------------|------|
| TCP                | 100                   | <1*                  | -    |
| TCMX               | 20                    | 98                   | 1.8  |
| Hexachlorobenzene  | 5                     | 89                   | 1.3  |
| alpha-BHC          | 5                     | 100                  | 1.4  |
| gamma-BHC          | 5                     | 99                   | 1.3  |
| beta-BHC           | 5                     | 102                  | 3.6  |
| delta-BHC          | 5                     | 103                  | 4.6  |
| Heptachlor         | 5                     | 91                   | 4.1  |
| Aldrin             | 5                     | 98                   | 1.8  |
| Heptachlor epoxide | 5                     | 103                  | 1.1  |
| gamma-Chlordane    | 5                     | 105                  | 3.4  |
| alpha-Chlordane    | 5                     | 96                   | 2.0  |
| 4,4'-DDE           | 10                    | 99                   | 1.6  |
| Endosulfan I       | 5                     | 101                  | 2.3  |
| Dieldrin           | 10                    | 98                   | 2.6  |
| Endrin             | 10                    | 97                   | 5.1  |
| 4,4'-DDD           | 10                    | 97                   | 2.9  |
| Endosulfan II      | 10                    | 99                   | 3.9  |
| 4,4'-DDT           | 10                    | 97                   | 2.2  |
| Endrin aldehyde    | 10                    | 92                   | 4.6  |
| Methoxychlor       | 50                    | 95                   | 4.1  |
| Endosulfan sulfate | 10                    | 100                  | 2.1  |
| Endrin ketone      | 10                    | 98                   | 1.3  |
| DCB                | 20                    | 101                  | 2.4  |

\*Meets requirement of <5% breakthrough.

# EPA 8081 Pesticides Extract Cleanup



Overlay of eight samples for each lot

| Peaks                           | ts (min) |
|---------------------------------|----------|
| 1. 2,4,5,6-Tetrachloro-m-xylene | 5.65     |
| 2. Biphenyl, decachloro-        | 10.00    |

**Column** Rtx-CLPesticides, 30 m, 0.32 mm ID, 0.32 µm (cat.# 11141)  
**Standard/Sample** 2,4,5,6-Tetrachloro-meta-xylene (cat.# 32027)  
Decachlorobiphenyl (BZ #209) (cat.# 32029)

**Diluent:** Hexane  
**Conc.:** 5 ng/µL

**Injection**  
**Inj. Vol.:** 4 µL pulsed splitless  
**Liner:** Topaz single taper inlet liner w/wool (cat.# 23303)  
**Inj. Temp.:** 250 °C  
**Pulse Pressure:** 35 psi (241.3kPa)  
**Pulse Time:** 0.74 min  
**Purge Flow:** 50 mL/min

**Oven**  
**Oven Temp.:** 70 °C (hold 0.5 min) to 320 °C at 25 °C/min (hold 2 min)

**Carrier Gas** He, constant flow  
**Flow Rate:** 3.5 mL/min

**Detector** Micro-ECD @ 330 °C  
**Make-up Gas Flow Rate:** 60 mL/min

**Make-up Gas Type:** N<sub>2</sub>  
**Data Rate:** 50 Hz

**Instrument** Agilent 7890B GC

**Sample Preparation** Conditioned a Resprep FL + CarboPrep Plus SPE cartridge (cat.# 28899) by adding one cartridge volume of hexane:acetone (90:10) and letting it stand for five minutes before drawing the solvent down to frit level. Then, 1 mL of sample extract was loaded onto the cartridge and eluted with 9 mL of hexane:acetone (90:10).



# Conclusion

The combined effects of both Florisil and GCB in a serial format provide remarkably increased removal of performance deteriorating coextracted matrices while still retaining the attributes and solvent elution ability of a standard Florisil cartridge thereby enabling better quality data and greater profitability

# Questions?

