

Isotopic “Fingerprinting” of PFAS – via the Thermo Scientific™ Orbitrap Exploris™ Isotope Solutions

Taylor Graham, Brett Davidheiser, Nils Kuhlbusch, Andreas Hilkert, Dieter Juchelka

National Environmental Monitoring Conference 2025

August 7, 2025

 The world leader in serving science



Outline

- Thermo Fisher's magnetic sector portfolio
- Why use an orbitrap for natural abundance isotope measurements?
- Proof of concept: natural abundance ratios on an orbitrap
- Measurement technique
- PFAS data
- Questions



The current Thermo Scientific™ gas IRMS portfolio

We currently offer a scaled portfolio of MS and sample introduction solutions



- **Thermo Scientific™ DELTA Q™ IRMS**
- Sensitivity up to 800 M/I
- Mass range up to m/z 96
- Accommodates up to 10 collectors



- **Thermo Scientific™ 253 Plus™ 10 kV IRMS**
- High-precision isotope analysis
- Long term stability and robustness
- Complete automation for ease-of-use



- **Thermo Scientific™ Ultra™ HR-IRMS**
- Double-focusing magnetic-sector mass analyzer
- Variable collector array
- Dual viscous flow inlet systems
- Optional (Multi) Ion Counting and RPQ

Oxyanions: refine understanding of element cycling


Atmosphere



Soil



Water





Oxyanions are important chemical constituents of virtually every environment on Earth.

What are current isotope ratio techniques missing?

A wish list

- Direct measurement of intact molecules
- Access to liquid, polar samples
- Higher sensitivity
- Less sample preparation, eliminate solid $\rightarrow$ gas process
- Position specific isotope analysis




Nitrogen **2021**, 2, 196–217. <https://doi.org/10.3390/nitrogen2020013>

Water **2020**, 12(1), 138; <https://doi.org/10.3390/w12010138>

An opportunity to bridge a gap exists


Combine technology developed for Life Science with approaches from stable isotope geochemistry




Mass spectrometers for isotope ratio analysis

Classic isotope-ratio MS N_2 , O_2 , N_2O


- **Indirect:** Conversion into **simple gases**
- Combining isotopologues on **few signals**



Thermo Scientific Orbitrap™ MS



- **Direct:** Intact isotopologues
- Separation of **all species** by HRMS



Orbitrap Exploris Isotope Solutions




Thermo Scientific™ Orbitrap Exploris™ Isotope Solutions includes:

Thermo Scientific™
Orbitrap Exploris™
120/240/480 MS



Dual Syringe Inlet

Data evaluation Package
for Isotope Ratio MS




Optional Thermo
Scientific™ Vanquish™
Neo UHPLC System



What's inside the Orbitrap Exploris MS

Essential components for isotope ratio analysis




Important parameters for Orbital analyzer



Isotope ratio analysis


- Resolution:
 - Orbitrap analyzes a specific number of ions (ion package) per scan
 - Every scan results in a mass spectrum
 - Higher resolution requires longer scan time:

	Resolution at m/z 200	Scans per second
NO ₃ ⁻	15,000	22
	30,000	12
HSO ₄ ⁻	45,000	10
	60,000	7
Organics	120,000	3
	240,000	1.5



Isotope ratio methodology for Orbitrap MS


- **Precision: maximizing ion counts**
 - Longer analysis times achieve greater precision.
 - Lower resolution setting requires less time for scanning.
 - Optimize mass range maximize the number of target ions.
- **Accuracy: Sample/Standard comparison**
 - NO_3^- Analysis of signals (reference/sample) for 6–9 mins



Dual Syringe Inlet

Sample introduction technique

- **Dual Syringe Inlet System** utilizing the diverter valve option of the Orbitrap Exploris MS
- Direct infusion of sample (**50 μM** in MeOH) with a flow rate of **4 $\mu\text{L}/\text{min}$** via a syringe pump; sample reference comparison by switching of a diverter valve.




- **In-flow Injection** - Vanquish Neo UHPLC System
- Loop-injection of **20-30 µL** sample (**50 µM** in MeOH) by the Autosampler into a flow of **4 µL/min** of MeOH.



Automated In-Flow Injection with Syringe Reference


Sample introduction technique

- **In-flow Injection** - Vanquish Neo UHPLC System
- Syringe injected reference material to control for ESI instrument drift. Increases orbitrap uptime while loop is washed out.




Orbitrap for isotopes: workflow for nitrate

Nitrate simplest model (4 atoms with isotope species: ^{14}N , ^{15}N , ^{16}O , ^{17}O , ^{18}O)



ESI-Orbitrap for isotopes – methodology

‘M0’ experiment




- M0 peak as basepeak

	<i>m/z</i>	Isotopologue	Abundance	Percentage
M0	61.9884	¹⁴ N ¹⁶ O ₃	989242	98.9
M+1	62.9854	¹⁵ N ¹⁶ O ₃	3637	0.36
	62.9926	¹⁴ N ¹⁷ O ¹⁶ O ₂	1127	0.11
M+2	63.9896	¹⁵ N ¹⁷ O ¹⁶ O ₂	4.1	
	63.9926	¹⁴ N ¹⁸ O ¹⁶ O ₂	5951	0.60
	63.9968	¹⁴ N ¹⁷ O ₂ ¹⁶ O	0.4	
M+3	64.9897	¹⁴ N ¹⁸ O ¹⁶ O ₂	21.9	
	64.9938	¹⁵ N ¹⁷ O ₂ ¹⁶ O	< 0.1	
	64.9968	¹⁵ N ¹⁷ O ¹⁸ O ¹⁶ O	4.5	
	65.0010	¹⁴ N ¹⁷ O ₃	< 0.1	
M+4	65.9939	¹⁵ N ¹⁸ O ¹⁶ O ₂	< 0.1	
	65.9969	¹⁵ N ¹⁸ O ₂ ¹⁶ O	11.9	
	65.9981	¹⁵ N ¹⁷ O ₃	< 0.1	
	66.0011	¹⁵ N ¹⁷ O ₂ ¹⁸ O	< 0.1	

ESI-Orbitrap for isotopes – methodology

‘noM0’ experiment



- ¹⁵N or ¹⁸O peak as base peak

	<i>m/z</i>	Isotopologue	Abundance	Percentage
M0	61.9884	¹⁴ N ¹⁶ O ₃	989242	98.9
M+1	62.9854	¹⁵N ¹⁶ O ₃	3637	33.81
	62.9926	¹⁴ N ¹⁷O ¹⁶ O ₂	1127	10.48
M+2	63.9896	¹⁵ N ¹⁷ O ¹⁶ O ₂	4.1	
	63.9926	¹⁴ N ¹⁸O ¹⁶ O ₂	5951	55.32
	63.9968	¹⁴ N ¹⁷ O ₂ ¹⁶ O	0.4	
M+3	64.9897	¹⁵N¹⁸O ¹⁶ O ₂	21.9	0.20
	64.9938	¹⁵ N ¹⁷ O ₂ ¹⁶ O	< 0.1	
	64.9968	¹⁴ N ¹⁷O¹⁸O ¹⁶ O	4.5	0.04
	65.0010	¹⁴ N ¹⁷ O ₃	< 0.1	
M+4	65.9939	¹⁵ N ¹⁸ O ¹⁶ O ₂	< 0.1	
	65.9969	¹⁴ N ¹⁸O₂ ¹⁶ O	11.9	0.11
	65.9981	¹⁵ N ¹⁷ O ₃	< 0.1	
	66.0011	¹⁴ N ¹⁷ O ₂ ¹⁸ O	< 0.1	

Dual Inlet measurements of nitrate




‘M0’ experiment:

- Three nitrate reference materials available from USGS:[1]

	$\delta^{15}\text{N}_{\text{AIR}}$	$\delta^{18}\text{O}_{\text{VSMOW}}$	$^{17}\text{O}_{\text{VSMOW}}$
USGS32	+ 180 ‰	+ 25.7 ‰	-
USGS34	- 1.8 ‰	- 27.9 ‰	- 14.8 ‰
USGS35	+2.7 ‰	+ 57.5 ‰	+ 51.5 ‰

Quality control: $^{15}\text{N}/^{14}\text{N}$ isotope ratio analysis

Block	Description	Ratio (¹⁵ N/ ¹⁴ N)	δ ¹⁵ N _{USGS32/Air} [‰]	δ ¹⁵ N _{USGS32/Air} [‰]	Std. Dev.		
1	Reference (USGS35)	0.00430	179.5	179.6	0.4		
2	Sample (USGS32)	0.00506					
3	Reference (USGS35)	0.00430					
4	Sample (USGS32)	0.00506	179.3				
5	Reference (USGS35)	0.00430	180.1				
6	Sample (USGS32)	0.00506					
7	Reference (USGS35)	0.00430					



[1] Böhlke, J.K. ; et al. Rapid Commun. Mass Spectrom., 2003, 17, p. 1835–1846.


Referencing scheme for nitrate

Analysis of international standards for nitrate

- Isotope ratio data of 'M0' and 'noM0' experiments using N11 as a working standard:[1]

	$\delta^{15}\text{N}_{\text{AIR}}$	$\delta^{18}\text{O}_{\text{VSMOW}}$	$\delta^{17}\text{O}_{\text{VSMOW}}$
USGS32	+ 180 ‰	+ 25.7 ‰	-
USGS34	- 1.8 ‰	- 27.9 ‰	- 14.8 ‰
USGS35	+2.7 ‰	+ 57.5 ‰	+ 51.5 ‰


Isotope ratio data of six of nitrate's most abundant isotopologues using the In-Flow Injection setup.



[1] Böhlke, J.K. ; et al. Rapid Commun. Mass Spectrom., 2003, 17, p. 1835–1846.

Nitrate analysis of environmental samples

'M0' experiment




Dual Syringe Inlet methodology

- Sample prep: “Dilute and Shoot”
 - Dilution to **1 μM (62 pg $\text{NO}_3^-/\mu\text{L}$)** for high salinity samples
- Results $\delta^{15}\text{N}$ and $\delta^{18}\text{O}$
 - Precision: < 1‰
 - Difference Orbitrap – Expected: < 1‰
- Takeaway
 - ESI tolerates up to **100-fold Cl⁻** load
 - Ground waters diluted up to **1/500**

Carbon stable isotope analysis of PFAS sources


- Carbon isotope analysis of PFAS, PFOA and FTOH were acquired from multiple vendors and where possible, from multiple lots.
- There was variation in the carbon isotope values among different compounds, manufacturers, and batches.
 - Ex. PFOA Syn $\delta^{13}\text{C}$ lots differed ($-28.9 \pm 0.3\text{‰}$, and $-46.1 \pm 0.7\text{‰}$).
- Variability among mean carbon isotope values makes EA-IRMS a potentially powerful tracking tool for PFAS contaminants.




- Dombrowski et al (2025). *Environmental Science & Technology Letters* <https://doi.org/10.1021/acs.estlett.5c00021>


FC(F)(F)C(F)(F)C(F)(F)C(F)(F)C(F)(F)C(F)(F)C(F)(F)C(F)(F)C(F)(F)C(=O)O

8 μ M solution of PFOA run on Orbitrap




Orbitrap Method Validation






PFOA

Chemical Formula: $C_8HF_{15}O_2$
Exact Mass: 413.97370




PFOS

1763-23-1
Chemical Formula: $C_8HF_{17}O_3S$
Exact Mass: 499.93749




PFBA

Chemical Formula: $C_4HF_7O_2$
Exact Mass: 213.98648




TFA

Chemical Formula: $C_2HF_3O_2$
Exact Mass: 113.99286




PFMOAA 674-13-5

Chemical Formula: $C_3HF_5O_3$
Exact Mass: 179.98458




HFPO-DA
13252-13-6

Chemical Formula: $C_3F_7O^*$
Exact Mass: 184.98374



Chemical Formula: $C_3HF_4O_2^*$
Exact Mass: 144.99127

HFPO-DA in-source fragmentation very common




PFBS

375-73-5
Chemical Formula: $C_4HF_9O_3S$
Exact Mass: 299.95027


Conclusions

- PFAS compounds are possible to isotopically fingerprint using Orbitrap IRMS
- Preliminary results show excellent agreement between traditional EA-IRMS and novel Orbitrap IRMS measurements
- Workflows on Orbitrap IRMS with PFAS compounds in solution are less damaging to instrumentation and provide faster throughput
- Orbitrap IRMS can be a powerful tool in source identification of PFAS contaminated sites.



Demo Lab at CU Boulder

- Orbitrap 240 with Vanquish Neo UHPLC
- Thermo Fisher: Brett Davidheiser
- PIs: Sebastian Kopf and Cajetan Neubauer
- Post Doc: Kristyna Kantnerova (now Assistant Prof at UCT Prague)



Thank you to our collaborators



Cajetan Neubauer and Kristýna Kantnerová for close collaboration and help

Joel Savarino, John K. Böhlke and Shohei Hattori for providing sample and reference materials

Thermo Fisher Scientific :

Kyle Fort

Konstantin Aizikov

Alexander Makarov

Uwe Rickens

Tabiwang Arrey

Colin Wirth

Christian Klaas

Charles Cartwright

John K. Böhlke
Stan Mroczkowski



John Eiler
Nathan Dalleska



Sebastian Kopf



Xingchen Wang



Ralph Mead
Chad Lane
Paul Wojtal



Thank you

