LUMINULTRA

SARS-CoV-2 Quantification and Variants Detection in Wastewater Samples:

Lessons Learned from Analyzing Geographically Diverse Samples

Anda Quintero, Anusha Edupuganti, Douglas Gramajo, John Boswell, Julia Capenter, and Michael Waud LuminUltra Technologies

PPT Template- April 6, 2021

Fredericton Baltimore Florida About us

Founded in 1995, LuminUltra is a biological diagnostic testing company headquartered in Canada with operations in 6 countries. It is widely recognized globally as a leader in developing tests and reagents for environmental, industrial, and diagnostic monitoring and is a key supplier of COVID-19 clinical testing reagents to the Government of Canada. Customers in over 80 countries trust LuminUltra's technology, production reliability and history of customer service excellence to deliver their essential services in a safe state.

LuminUltra offices spread across the globe

3.5K

Customers from around the world leverage LuminUltra's testing solutions

Countries with customers that rely on LuminUltra's products and services

\$345M In customer value

delivered to date

25 Years of exceptional production reliability

and innovation

Clinical COVID-19 testing reagents per week being supplied to Canada

500K

Melbourne

LUMINULTRA

LuminUltra proudly serves some of the top companies in the world including:





Heerenveen

Paris

London

VEOLIA OUPONTE ECELAB







Challenges we are facing in SARS-CoV-2 wastewater monitoring since April 2020

- 1. Pandemic effects:
 - A. High volume of samples.
 - B. Short turn around time.
 - C. Unstable reagents and consumables supplies.
 - ---Pandemic effect put great pressure on lab personnel and quality management system.
- 2. Unknowns in evaluation of measurement uncertainty.
- 3. Market requirements of monitoring emerging variants with clinical relevance.
- 4. Educating clients about how to use their data.

Lessons Learned

- 1. Automation in the lab.
- 2. Properly assess sample processing spike-in recovery rates .
- 3. Constantly developing new assays to detect mutations and NGS solutions.
- 4. Knowledge sharing helps the public better understand their data and the importance of continuous monitoring effort.



Wastewater samples characters are consistent throughout the sampling period



SARS-CoV-2 Viral Concentration determined by two concentration methods

HA Filtration Concentration vs. Direct Extraction?





Theoretical LOD by dashed lines (3 copies/reaction, 95% detection rate, 100% recovery rate). For a typical 50ml sample, LOD is 1.5E3 copies/liter.

N1 (Copies/Liter) From Direct Extraction (DE) Samples



Theoretical LOD by dashed line (3 copies/reaction, 95% detection rate). For direct extraction of 1ml water sample, LOD is 7.2E4 copies/liter.

LUMINULTRA

9

DE samples results in higher N1 concentration than filter samples, but similar trend in N1 concentration throughout 3 sampling events.



LUMINULTRA

10

N1 and N2 results are highly correlated with each other in filter samples



N1 and N2 filter sample (copies/liter)







The Unknown: Workflow Recovery

Spike-in at sample processing vs. Spike-in at RNA extraction

12

LUMIN<mark>ULTRA</mark>

Workflow recovery is a black box





LUMIN<mark>ULTRA</mark>

Filter Samples: Positive correlation between sample processing spike-in (OC43) and extraction spike-in (MLR) recovery rates



- OC43 recovery is determined by comparison to a pre-determined value of genomic copies of OC43 in heat-inactivated viral fluid.
- MLR recovery rate is determined by Delta Ct method with reference to Ct of 0.5ng/ul MLR per PCR reaction.

DE samples: Lack of correlation between sample processing spike-in (OC43) and extraction spike-in (MLR) recovery rates



Extraction Spike-in Recovery

Quantification of Clinically Important SARS-CoV-2 Mutations by RT-ddPCR

Are wastewater samples suitable for variants monitoring by RT-ddPCR?



Percentage of DEL69/70 and E484K mutants in Composite Filter Samples



9



Percentage of DEL69/70 and E484K mutants in Composite Filter Samples

20



þ

Time_Point

- - - - - -

. . .

1e+03-

- - - - - - - -

.

_ _ _ _ _ _ _ .

- - - - - - - -

High N1 (Copies/Liter) from Site B and Site H Filter Samples



SARS-CoV-2 Variants Identification by Whole Genome Sequencing

NGS Workflow





Oxford Nanopore MinION Mk 1B device and sequencing workspace



Challenges of mutation identification in wastewater samples



Lessons Learned

- 1. Automation in the lab with magnetic beads extraction methods
 - \rightarrow dramatically increase sample throughput.
- 2. Sample processing spike-in recovery rate

 \rightarrow Workflow recovery measured by sample spike-in provide necessary information about method efficiency and sample nature.

- 3. Constantly developing new assays to detect mutations and NGS solutions. \rightarrow RT-ddPCR is a feasible and valid methodology for quantification of mutations.
- 4. Knowledge sharing helps the public better understand their data and the importance of continuous monitoring effort.

 \rightarrow Joint efforts from:

academic research, preprints publications, multi-lab comparison studies, monitoring networks, data sharing platforms, and local community participations.