



Apply US-EPA Wastewater & Ambient Water ATP Guideline to Validate Colitag™ for Enumeration of *E. coli* and Fecal Coliform Bacteria

Dr. Preetha Biswas, Director of Microbiology, R&D
Co-Author: Dr. Lei Zhang, Principal Scientist, Microbiology R&D

Outline

01

Background

02

Colitag - Overview

03

ATP Study Requirement

04

ATP Study Plan


05

Preliminary Data

Background

For Wastewater Testing


- Clean Water Act (CWA)
 - Established a national commitment to restore and maintain integrity of the nation's waters
 - Improving the health of rivers, lakes, and coastal waters
- National Pollutant Discharge Elimination System (NPDES) permit program
 - Section 402 CWA
 - Regulating point sources that discharge pollutants to waters



Background

For Wastewater Testing

- Microbiological water pollution
- Indicators of microbial water safety
- Commonly used indicator bacteria for water fecal contamination
 - Fecal coliform bacteria
 - *E. coli*
 - *Enterococci*





Colitag

Potable water testing



- Colitag™ originally designed for simultaneous detection of *E. coli* and Total Coliforms in drinking water.
- Colitag for Drinking Water testing is validated both for qualitative and quantitative. Works in presence/absence method, or most probable number-based methodology.

Colitag – Overview

Potable water testing

Colitag medium and method using multi-well MPNTray™ for enumeration of *E. coli* and fecal coliforms.

- Detection is based on the presence of two enzymes
 - β -galactosidase
 - Coliforms, yellow
 - β -glucuronidase
 - *E. coli*, fluorescence
- Quantification by most probable number (MPN) method using the number of positive versus negative wells
 - Multi-well tray (97-well)

ONPG - Coliforms (*B*-Galactosidase) → O - Nitrophenol

MUG - *E. coli* (*B*-Glucuronidase) → 4 - Methyl - Umbelliferone

Objective

For Wastewater Testing




- EPA Microbiological ATP Protocol and Methods to be followed as guideline to evaluate the performance of Colitag medium using multi-well MPNTray, for enumeration of *E. coli* and Fecal Coliform Bacteria in wastewater



ATP Study Requirement

Wastewater Testing

- Follow the Alternate Test Procedure (ATP) guidelines
- Collaboration with EPA coordinator and advisors for review and advice on study proposal – and approval of the method



ATP Study Requirement

Wastewater Testing

Side-by-Side Comparison

- Parallel testing of the ATP and the reference method
- Ten sewage samples from geographic diverse wastewater plant locations
 - Secondary effluent spiked into final sewage targeting to achieve a 20-200 MPN/100 mL sample
- Twenty replicates for each sample to be evaluated


Side-by-Side Comparison ATP Study Plan

Wastewater Testing

Sample source		Replicates	Target Analysis	Analysis method	Minimum comparability results
Type	Number				
Wastewater	10	20	FC	SM 9223	200
				Colitag™	200
	10	20	<i>E. coli</i>	SM 9223	200
				Colitag™	200


Side-by-Side Comparison Study Flow Chart

Wastewater Testing




Confirmation Flow Chart for Positive Results

Wastewater Testing



Confirmation Flow Chart for Negative Results

Wastewater Testing



Internal Preliminary Test – Recovery of *E. coli*

Colitag for Wastewater Testing

	Test 1, n=20		Test 2, n=20	
	Average counts per 100 mL	RSD (%)	Average counts per 100 mL	RSD (%)
Colitag	48	20.51	22	23.44
SM 9223B	55	14.23	26	26.66
1603	46	16.92	19	33.76

Verification	Colitag	SM 9223B	1603
Sensitivity	94.7%	86.4%	90%
Specificity	90.5%	94.4%	100%

Internal Preliminary Test – Recovery of FC

Colitag for Wastewater Testing

	Test 1, n=20		Test 2, n=20	
	Average counts per 100 mL	RSD (%)	Average counts per 100 mL	RSD (%)
Colitag	23	21.00	78	18.99
SM 9223B	13	33.56	75	17.28
SM 9222D	24	34.012	66	14.30

Verification	Colitag	SM 9223B	SM 9222D
Sensitivity	100%	95.2%	100%
Specificity	100%	100%	100%

Internal Data: Inclusivity Testing

Colitag Specificity

Inclusivity Organisms		Inoculum (CFU/100 mL)	ONPG (18-24h, 35°C)	MUG (18-24h, 35°C)
• <i>Citrobacter freundii</i>	ATCC 6879	21	Positive	Negative
• <i>Citrobacter freundii</i>	ATCC 8090	9	Positive	Negative
• <i>Cronobacter sakazakii</i>	ATCC 29004	52	Positive	Negative
• <i>Enterobacter amnigenus</i>	ATCC 51818	43	Positive	Negative
• <i>Escherichia coli</i>	ATCC 8739	18	Positive	Positive
• <i>Escherichia coli</i>	ATCC 25922	26	Positive	Positive
• <i>Escherichia coli</i>	CDC 984	88	Positive	Positive
• <i>Escherichia coli</i>	ATCC 11775	50	Positive	Positive
• <i>Escherichia coli</i>	ATCC 11229	20	Positive	Positive
• <i>Hafnia alvei</i>	ATCC 51815	42	Positive	Negative
• <i>Hafnia paralvei</i>	FSL-W4-0268	47	Positive	Negative
• <i>Klebsiella oxytoca</i>	ATCC 2170	28	Positive	Negative
• <i>Klebsiella pneumoniae</i>	ATCC 13883	21	Positive	Negative
• <i>Klebsiella pneumoniae</i>	ATCC 27736	16	Positive	Negative
• <i>Serratia marcescens</i>	ATCC 8100	20	Positive	Negative

Internal Data: Exclusivity Testing

Colitag Specificity

Exclusivity Organisms		Inoculum (CFU/100 mL)	ONPG (48h, 35°C)	MUG (48h, 35°C)
• <i>Aeromonas hydrophilia</i>	ATCC 35654	9.1 x 10 ⁵	Negative	Negative
• <i>Aeromonas hydrophilia</i>	ATCC 7966	4.0 x 10 ⁶	Negative	Negative
• <i>Acinetobacter baumannii</i>	ATCC 19606	5.4 x 10 ⁵	Negative	Negative
• <i>Alcaligenes faecalis</i>	ATCC 925	1.0 x 10 ⁶	Negative	Negative
• <i>Enterococcus faecalis</i>	ATCC 29212	1.9 x 10 ⁶	Negative	Negative
• <i>Kocuria rizophilia</i>	ATCC 9341	6.7 x 10 ⁵	Negative	Negative
• <i>Pseudomonas aeruginosa</i>	ATCC 9027	2.3 x 10 ⁶	Negative	Negative
• <i>Pseudomonas aeruginosa</i>	ATCC 10145	2.8 x 10 ⁶	Negative	Negative
• <i>Staphylococcus aureus</i>	ATCC 25923	3.3 x 10 ⁵	Negative	Negative
• <i>Bacillus subtilis</i>	ATCC 9372	2.7 x 10 ⁴	Negative	Negative
• <i>Bacillus subtilis</i>	ATCC 6633	4.3 x10 ⁵	Negative	Negative

Side-by-Side Comparison ATP Study Plan - Recap

Wastewater Testing

Sample source		Replicates	Target Analysis	Analysis method	Minimum comparability results
Type	Number				
Wastewater	10	20	FC	SM 9223	200
				Colitag™	200
	10	20	<i>E. coli</i>	SM 9223	200
				Colitag™	200

ATP Study Side-by-Side Comparison

Wastewater Testing

Review of study result:

- Mean recovery for each matrix

$$\% \text{ Recovery} = \frac{\text{Result} - \text{Background}}{\text{Spike}} * 100 \%$$

Where,

Result = the amount recovered from the sample after spiking,

Spike = the estimated amount spiked into the sample, and

Background = the estimated background amount measured in the sample prior to spiking.

- Precision: RSD

$$RSD = \frac{SD}{\text{Mean}} * 100 \%$$

Where,

SD = the standard deviation of all recovered amounts

Mean = the mean of all recovered amounts

ATP Study Side-by-Side Comparison

Wastewater Testing

Review of study result:

- False negative rates

		Method		Total
		New	Reference	
Result	True +	TP ₁	TP ₂	TP ₁ + TP ₂
	False -	FN ₁	FN ₂	FN ₁ + FN ₂
Total		TP ₁ +FN ₁	TP ₂ + FN ₂	TP ₁ + TP ₂ + FN ₁ + FN ₂

- False positive rates

		Method		Total
		New	Reference	
Result	False +	FP ₁	FP ₂	FP ₁ + FP ₂
	True -	TN ₁	TN ₂	TN ₁ + TN ₂
Total		FP ₁ + TN ₁	FP ₂ + TN ₂	FP ₁ + FP ₂ + TN ₁ + TN ₂

- QC Acceptance Criteria-Based
Method Blank, Positive Control and Media Sterility

Thanks to the EPA Team

Thank you
for listening

Get in touch

www.NEOGEN.com

