



# **Measurement of Gross Alpha: Application, Interpretation and New Method Options**

Glynda Smith

US Environmental Protection Agency  
Office of Groundwater and Drinking Water  
Technical Support Center  
Cincinnati, Ohio




## Disclaimer

The information in this presentation has been reviewed and approved for public dissemination in accordance with U.S. Environmental Protection Agency (EPA). The views expressed in this presentation are those of the author(s) and do not necessarily represent the views or policies of the Agency. Any mention of trade names or commercial products does not constitute EPA endorsement or recommendation for use.

## Some fundamentals

Chemical bonding/breaking/ ionic interactions – electron shells




Radioactive decay – nucleus

© 2012 Encyclopædia Britannica, Inc.

## Some fundamentals

- All elements above atomic number 83 (bismuth) are radioactive
- Release of alpha particle (helium nucleus) results in a decay product with reduced atomic number (-2) and mass (-4).



Source: opentextbc.ca



## Some fundamentals

- Activity ~ degree of radioactivity = decays/time
  - picocurie, pCi = 2.22 disintegrations/min
  - becquerel, Bq = 27.03 pCi
- Half-life ( $t_{1/2}$ ) = time it takes the number of radioactive atoms to be reduced to one-half the original amount
  - Each radioactive isotope has its own characteristic half-life
- Activity is inversely proportional to half-life  
Shorter half-life means higher activity



## Background


- Many groundwaters contain measurable amounts of naturally occurring alpha particle-emitting radionuclides
- Gross alpha analysis provides a measurement estimate of alpha particle activity
- Gross alpha analysis cannot discriminate among all possible alpha emitters or their energies



# Background: Decay Chains

Most naturally occurring alpha emitters arise from two decay chains. Uranium-238 decay chain includes U-238, U-234, Th-230, Ra-226, Rn-222, Po-218, Po-214, Po-210

## The Uranium-238 Decay Chain




Source: U.S. Geological Survey (USGS)



## Background: Decay Chains

The Thorium-232 decay chain includes Th-232, (Ra-228 a beta-emitter), Th-228, Ra-224, Rn-220, Po-216, Bi-212, Po-212, Pb-208, and Tl-208.

### The Thorium-232 Decay Chain



Source: U.S. Geological Survey (USGS)






## Assumptions

- Cannot calibrate detector for every possible alpha particle emitter
  - Calibrate with a “reference radionuclide”
  - Approved drinking water reference radionuclides (40 CFR 141.25(a))
    - Gross alpha evaporation methods (e.g. EPA Method 900.0): Th-230, U-nat (not Am-241)
    - Gross alpha co-precipitation methods (e.g. Standard Method 7110 C): Th-230, U-nat, and Am-241
- Gross alpha analysis implicitly assumes that all alpha particle emitters are detected with the same efficiency as the reference radionuclide.



### Example Alpha Energies, MeV




MeV: Mega-electron volts



## Limitations - Geometry

- Geometry is critical!



Source: Arndt and West. 2007. *Health Physics*, 92, 148-156



## Limitations - Geometry

- Consistent preparation of uniform sample residues is key!
  - Alpha particle energy emission reduced by any kind of barrier
  - Good sensitivity only possible at low total solids levels
    - Limited to samples that yield  $< 100$  mg solid residue in a  $20\text{ cm}^2$  counting planchet (i.e.  $< 5\text{ mg/ cm}^2$ )
  - “Patchy” residues inhibit alpha particle energy emission more than uniform residues
  - Emissions from center of planchet are more likely to be captured by detector than those from edges



## Gross Alpha Variability – Timing Events

When is the best time to count? Depends on what you want to know!

- Radium-224 (half-life 3.66 days)
  - Samples containing radium-228 activity can contain comparable radium-224 activity at the time of collection (remember the thorium-232 decay chain – radium-224 is one of the progeny of radium-228)
  - Radium-224 produces four short-lived alpha particle emitters: Rn-220, Po-216, Bi-220, and Po-212
  - Radium-224 activity is rarely determined (except in some states), but if a sample is prepared and counted within one week of collection, radium-224 and its progeny can have a significant impact on gross alpha activity
  - If radium-224 contribution to gross alpha activity is not required by your state, hold samples 2-3 weeks before preparing/counting to allow radium-224 and its progeny to decay away.



## Gross Alpha Variability – Timing Events


40 CFR 141.66(c): Gross alpha particle activity MCL is 15 pCi/L (including radium-226, but excluding radon and uranium)

- Screen water samples for radium-226 at levels that may exceed the combined radium MCL of 5 pCi/L
- The most accurate screening for radium-226: delay sample preparation for 2-3 weeks (allow radium-224 and progeny to decay), then prepare and count *right away*.
- If radium-226 is present in a sample, the time between sample preparation and counting can significantly impact gross alpha. Why? Radon!
- Radon is removed when the sample is prepared. But sample residue forms a thin film. As radium-226 decays, radon-222 is formed and trapped in the solid residue.



## Gross Alpha Variability – Timing Events

As radon-222 and its short-lived alpha-emitting progeny polonium-218 and polonium-214 build up within the prepared evaporated sample, alpha activity increases. In 3-4 weeks after sample preparation, activity quadruples!



Eaton and Shannon. 2015. AWWA. Radionuclide Rule Compliance: Utility Guidance on Analytical Methods



## Gross Alpha Variability – Timing Events

EPA Method 900.0, Sect. 8.4: After drying (the sample) “alpha counting should be delayed at least 72 hours until equilibrium has occurred.” No consensus on what ‘equilibrium’ means.

If radium-226 is present (and radium-224 is not considered), 72 hours allows two things to happen:

- When sample is prepared, radon-222 is removed but alpha-emitting progeny remain. Decay of the unsupported progeny will occur within ~ 8 hours
- Partial ingrowth of radon-222 and its supported progeny from radium-226 begins

Specifying that prepared samples be held at least 72 hours can result in inter-laboratory inconsistencies. Consider what happens if one lab counts after 72 hours and another lab waits a week or two.





## What's with Gross Alpha and Uranium?

- Uranium is regulated as a heavy metal (MCL = 30  $\mu\text{g/L}$ ), but it is also an alpha particle-emitting radionuclide
- Radionuclide Rule specifies conversion factor of 0.67 pCi/ $\mu\text{g}$  to convert uranium mass concentration in  $\mu\text{g/L}$  to an estimated activity concentration in pCi/L
  - 0.67 pCi/ $\mu\text{g}$  based on activity ratio in uranium-bearing minerals (not water)
  - Yields a conservative estimation of uranium activity from a public health protection standpoint, i.e. calculated uranium alpha emission activity will be biased low, resulting in a higher adjusted gross alpha (gross alpha minus uranium contribution)
  - If adjusted gross alpha based on converting uranium mass to activity exceeds 15 pCi/L MCL, determine a more accurate uranium activity using a radiochemical method (e.g. alpha spectrometry)



## What's Being Done to Reduce Gross Alpha Variability?

- Still a screening method and some limitations will exist (e.g. calibration reference radionuclide)
- But variability in some factors can be controlled:
  - Revising the 900-Series EPA Methods
  - Approving new methodologies for gross alpha analysis



## EPA Method 900.0, Rev. 1.0 (Feb. 2018)

**EPA Doc. 815-B-18-002;  
Approved Oct. 12, 2018 (83 FR  
51636)**

- Updated for simultaneous alpha/beta counting and addresses spillover (crosstalk)
- Details efficiency calibration and QC specifications
- Details the importance of timing events and removes the minimum 72-hour hold time between sample preparation and counting



Method 900.0, Revision 1.0: Gross Alpha and Gross Beta  
Radioactivity in Drinking Water



## **ASTM D7283-17 and Standard Method 7110 D-17**

**Approved July 27, 2017 (82 FR 34861)**

- Gross alpha and gross beta by liquid scintillation counting (LSC)
  - Simultaneous alpha/beta counting with pulse height discrimination
  - Consensus organizations followed the EPA Drinking Water Radiochemistry Alternate Test Procedure (ATP) Protocol to set up multi-lab validation study, QC and acceptance criteria
  - Newer scintillation cocktails provide homogeneous sample dissolution
    - Minimizes matrix solids effects
    - LSC provides higher detection efficiencies



## SDWA Radionuclide Detection Limit

EPA Doc. 815-B-17-003; Published April 2017

- DL defined at 40 CFR 141.25(c)
- SDWA DL  $\neq$  MDA, MDC or MDL
- Derives equation for reporting:

$$DL = \frac{\frac{1.96^2}{2t_s} \times \left[ 1 + \sqrt{1 + \frac{4t_s^2}{1.96^2} \times B \times \left( \frac{1}{t_s} + \frac{1}{t_B} \right)} \right]}{\varepsilon \times V \times 2.22}$$



## Summary

- Several sources of variability can bias gross alpha results:
  - Reference radionuclide used for calibration
  - Different energies of alpha emitters
  - Uniformity of prepared sample residues
  - Timing intervals between collection and preparation of sample, and between preparation and counting of sample
- Improving consistency through new methods:
  - Current analytical technology
  - Quality control criteria
  - Clearly discuss the importance of timing events



## References for Help

- Approved analytical drinking water methods  
<https://www.epa.gov/dwanalyticalmethods/approved-drinking-water-analytical-methods>.
  - Select Radionuclides
  - Table of approved methods with links to the approved versions (including EPA Method 900.0, Rev. 1.0). Table is current as of April 2019
- Radionuclide DL document and radiochemistry training modules  
<https://www.epa.gov/dwlabcert>.
  - Under “Certification, Monitoring and Analysis Requirements” select the Procedure for Safe Drinking Water Act Program Detection Limits for Radionuclides
  - Under “Certification Process and Training” select Radiochemistry Training for Drinking Water Compliance



# Questions

