



An Inter-laboratory Study on EPA Methods 537.1 and 533 for Potable and Non-Potable Water PFAS Analyses

Joshua Whitaker, Marnellie Ramos, and Yongtao Li, Eurofins Eaton Analytical, LLC


Robert B. Hrabak, Eurofins TestAmerica

Christine Ratcliff and Charles Neslund, Eurofins Lancaster Laboratories Environmental, LLC



More than 3000 per- and polyfluoroalkyl substances (PFASs) are, or have been, on the global market.....


Courtesy to Wang, DeWitt, Higgins & Cousins, *Environ. Sci. Technol.* 2017, 51(5), 2508-2518.



* PFASs in RED are those that have been restricted under national/regional/global regulatory or voluntary frameworks, with or without specific exemptions (for details, see OECD (2015), Risk reduction approaches for PFASs. <http://oe.cd/iAN>).
 ** The numbers of articles (related to all aspects of research) were retrieved from SciFinder® on Nov. 1, 2016.

Figure 1. "Family tree" of PFASs, including examples of individual PFASs and the number of peer-reviewed articles on them since 2002 (most of the studies focused on long-chain PFCAs, PFSAs and their major precursors.).

Types of PFAS




Current PFAS Method Summary



Method	EPA 537/537.1	EPA 533	SW-846 EPA 8328	SW-846 EPA 8327	ASTM D7979-17	DoD/DoE QSM 5.3	Lab Methods	TOP Assay	TOF Assay
Techniques	RP-SPE LC/MS/MS IS CAL	WAX-SPE LC/MS/MS ID CAL	WAX-SPE LC/MS/MS ID CAL	1:1 MeOH LC/MS/MS Ext CAL	1:1 MeOH LC/MS/MS Ext CAL	WAX-SPE LC/MS/MS ID CAL	WAX-SPE LC/MS/MS ID CAL	Oxidation SPE LC/MS/MS	Combustion IC
Availability	2009 (v1.1)/ 2020 (v2)	2019	Draft	2019	2017	2019	Varies	Varies	Varies
Type of Water	Potable	Potable	Non-potable	Non-potable	Non-potable	Non-potable	Non-potable	Non-potable	Non-potable
Type of PFAS	14/18 analytes	25 analytes	24 analytes?	24 analytes	14+7 analytes	25 analytes	~40 analytes	All perfluoro- alkyl acids	All organo- fluorine
Type of analysis	Quantitative	Quantitative	Quantitative	Quantitative	Quantitative	Quantitative	Quantitative	Semi-quant. Screening	Semi-quant. Screening
Pros & Cons	Validated Limited S.C.	Validated Limited L.C.	Draft	Validated Faster Higher RLs	Validated Faster Higher RLs	Validated	More PFAS Inconsistent	All perfluoro- alkyl acids	Surrogate of all fluorine compounds
MRLs	~2 ng/L	~2 ng/L	~2 ng/L	10-50 ng/L	10-50 ng/L	~2 ng/L	~2 ng/L	ng-µg/L	µg-mg F/L

PFAS Have Been around for a Long Time.



Drinking Water PFAS Regulations



➤ USEPA

- 2013-2015 UCMR3 for 6 PFAS with MRLs of 10-90 ng/L by EPA Method 537
- Health Advisories (2016): PFOA/PFOS or PFOA + PFOS = 70 ng/L
- 2023-2025 UCMR5 for 29 PFAS with MRLs of 2-8 ng/L by EPA Methods 537.1 and 533
- MCLs for PFOA, PFOS, and other PFAS?

➤ States with DW Regulations as June 30, 2021

- Established / interim / proposed MCL: NJ, NH, VT, MA, NY, WI, ME
- Established NL: CA
- Proposed HBV / AL / TL: MI, MN, OH, RI, NC, OR, CO, IL

➤ Examples

- Individual PFAS NJ MCL: PFNA = 13 ng/L, PFOA = 14 ng/L, PFOS = 13 ng/L
- Sum of PFAS MA MCL: PFOA + PFHxS + PFOS + PFHpA + PFNA + PFDA = 20 ng/L

Proposed UCMR5 Timeline of Activities




2022	2023	2024	2025	2026
<p>Pre-sampling Activity by EPA</p> <ul style="list-style-type: none"> • Manage Lab Approval Program • Organize Partnership Agreements and State Monitoring Plans • Begin PWS SDWARS registration/inventory • Review GWRMP submittal • Conduct outreach/trainings 	<p>Sampling Period</p> <p>EPA Implementation Activities</p> <ul style="list-style-type: none"> • Provide compliance assistance • Implement small system monitoring • Post data quarterly to NCOD <p>PWS Sample Collection; Laboratory Analysis; Reporting</p> <ul style="list-style-type: none"> • All large systems serving more than 10,000 people; • All small systems serving between 3,300 and 10,000 people; • 800 small systems serving fewer than 3,300 people 			<p>Post-sampling Activity</p> <p>PWSs, Laboratories</p> <ul style="list-style-type: none"> • Complete resampling, as needed • Conclude data reporting <p>EPA</p> <ul style="list-style-type: none"> • Complete upload of UCMR 5 data to NCOD

Proposed UCMR5 29 PFAS (2023-2025)



EPA 533 (2-5 ng/L)	PFBA	PFPeA	PFHxA	PFHpA UCMR3, 10 ng/L	PFOA UCMR3, 20 ng/L
	PFNA UCMR3, 20 ng/L	PFDA	PFUnA	PPDoA	PFBS UCMR3, 90 ng/L
	PFPeS	PFHxS UCMR3, 30 ng/L	PFHpS	PFOS UCMR3, 40 ng/L	PFEESA
	4:2 FTS	6:2 FTS	8:2 FTS	HFPO-DA	ADONA
	9Cl- PF3ONS	11Cl- PF3OUdS	PFMBA	PFMPA	NFDHA
EPA 537.1 (5-8 ng/L)	MEtFOSAA	NMeFOSAA	PFTTrDA	PFTeDA	

EPA 537.1 vs. EPA 533



Method	EPA 537.1	EPA 533
250 mL samples	14 days Trizma pH 6.5 - 7.5	28 days ammonium acetate pH 6 - 8
1 mL extracts	28 days 96% MeOH/water,	28 days 80% MeOH/water
IS / IPS	Internal standards unextracted	Isotope performance standards unextracted
SS / IDA	Surrogate standards extracted	Isotope dilution analogues extracted
Calibration	Internal standard calibration	Isotope dilution calibration
QC	Typical	Typical

EPA 537.1 and 533 Applications



- **Drinking water utilities are often involved in two types of analyses.**
 - **Finished drinking water quality compliance**
 - **Source water assessments**

- **EPA Methods 537.1 (evolving from 537) and 533 were developed and validated as drinking water methods.**
 - **What are the main challenges for drinking water analyses?**
 - Relatively low EPA 537.1 surrogate NEtFOSAA-d5 recoveries
 - Field sample (FS) and field reagent blank (FRB) bottles switched

 - **Are these methods applicable for source water assessments?**
 - What are the main challenges for non-potable water analyses?
 - How to resolve these challenges?

Inter-Laboratory Study Water Matrices



RW1	Native (n = 1)	2 ng/L (n = 4)	RW2	Native (n = 1)	50 ng/L (n = 4)
DW1	Native (n = 2)	10 ng/L (n = 4)	DW2	Native (n = 2)	10 ng/L (n = 4)
GW1	Native (n = 2)	10 ng/L (n = 4)	GW2	Native (n = 2)	10 ng/L (n = 4)
SW1	Native (n = 2)	10 ng/L (n = 4)	SW2	Native (n = 2)	10 ng/L (n = 4)
WW1	Native (n = 2)	50 ng/L (n = 4)	WW2	Native (n = 2)	50 ng/L (n = 4)

Inter-Laboratory Study Water Matrices (Cont'd)



Matrix	DW1	DW2	GW1	GW2	SW1	SW2	WW1	WW2
pH	6.8	7.2	7.3	7.4	7.3	6.7	7.3	7.2
Free chlorine (mg/L)	0.86	0.86	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05
TOC (mg/L)	0.52	0.52	0.86	0.43	3.1	5.8	2.4	3.7
Total hardness as CaCO ₃ (mg/L)	401	563	257	256	278	122	412	324
Total alkalinity as CaCO ₃ (mg/L)	276	275	208	320	214	109	276	213
Chloride (mg/L)	147	289	38	37	33	13	239	182
Sulfate (mg/L)	60	60	28	27	39	4.6	94	54
Nitrate as nitrogen (mg/L)	< 1.0	< 1.0	1.7	2.1	1.8	< 1.0	11	13
HPC (MPN/mL)	NA	NA	311	372	650	1000	440	623

Mean Native Concentrations ± SD (ng/L)



Analyte	Method	DW1	DW2	GW1	GW2	SW1	SW2	WW1	WW2
PFBA	533	3.2 ± 0.3	3.2 ± 0.3	4.5 ± 0.3	4.6 ± 0.5	2.0 ± 0.3	2.8 ± 0.5	6.9 ± 1.8	7.8 ± 2.9
PFPeA	533	1.7 ± 0.1	1.7 ± 0.2	0.9 ± 0.1	0.8 ± 0.1	1.4 ± 0.3	1.4 ± 0.3	20.3 ± 1.7	15.4 ± 1.2
PFBS	533	2.4 ± 0.1	2.4 ± 0.2	5.4 ± 0.4	5.2 ± 0.4	1.3 ± 0.7	0.9 ± 0.5	3.5 ± 0.3	4.0 ± 0.6
	537.1	2.5 ± 0.2	2.5 ± 0.2	5.5 ± 0.2	5.5 ± 0.4	1.4 ± 0.8	0.9 ± 0.5	3.5 ± 0.7	4.0 ± 1.5
PFHxA	533	2.0 ± 0.1	2.0 ± 0.2	1.0 ± 0.1	0.9 ± 0.1	0.9 ± 0.6	1.0 ± 0.6	14.9 ± 1.7	15.9 ± 1.7
	537.1	1.6 ± 0.9	2.1 ± 0.2	0.7 ± 0.4	0.8 ± 0.5	1.0 ± 0.6	1.1 ± 0.7	14.7 ± 0.8	15.5 ± 1.1
PFHpA	533	1.8 ± 0.2	1.9 ± 0.2	1.0 ± 0.1	0.8 ± 0.1	ND	1.1 ± 0.2	1.8 ± 0.2	1.5 ± 0.3
	537.1	1.5 ± 0.9	1.5 ± 0.9	ND	ND	ND	1.1 ± 0.6	2.0 ± 0.2	1.5 ± 0.3

Mean Native Concentrations ± SD (ng/L)




Analyte	Method	DW1	DW2	GW1	GW2	SW1	SW2	WW1	WW2
PFHxS	533	1.2 ± 0.1	1.2 ± 0.1	0.8 ± 0.1	0.9 ± 0.1	0.7 ± 0.4	ND	2.9 ± 0.4	5.5 ± 0.6
	537.1	ND	1.0 ± 0.6	0.8 ± 0.5	0.7 ± 0.5	0.9 ± 0.6	ND	3.2 ± 0.3	5.9 ± 0.7
PFOA	533	9.6 ± 0.8	9.4 ± 0.7	3.0 ± 0.3	2.8 ± 0.2	0.8 ± 0.5	1.9 ± 1.1	4.0 ± 0.5	3.8 ± 0.6
	537.1	9.7 ± 0.8	10.1 ± 0.6	2.8 ± 0.2	2.8 ± 0.2	1.0 ± 0.6	2.7 ± 0.2	4.6 ± 0.5	3.9 ± 0.3
PFOS	533	ND	ND	ND	ND	1.1 ± 0.7	1.0 ± 0.6	2.5 ± 0.3	9.6 ± 0.6
	537.1	ND	ND	0.7 ± 0.6	0.7 ± 0.6	1.2 ± 0.9	1.1 ± 0.8	3.1 ± 0.9	9.4 ± 1.1
PFNA	533	ND	ND	ND	ND	ND	0.7 ± 0.5	ND	ND
	537.1	ND	ND	ND	ND	ND	ND	0.7 ± 0.3	ND

EPA 533 IDA Recoveries




- All multi-lab IDA recoveries were within 50-200%.
- The SW and WW matrices resulted in slightly low recoveries for long-chain IDA.
 - PFDA-13C6
 - PFUnA-13C7
 - PFDoA-13C2
- The WW matrices significantly affected IPS-PFBA-13C3 and IDA-PFBA-13C4 peak areas.



EPA 533 IDA %RSD




- Single lab RSD were $\leq 20\%$.
- The multi-lab IDA RSD were generally $\leq 20\%$.
- The SW and WW matrices resulted in relatively higher RSD.
 - PFDoA-13C2 = 23% for SW1
 - 4:2 FTS-13C2 = 23% for SW2
 - PFPeA-13C5 = 24% for WW2



EPA 533 Analyte Recoveries




- All multi-lab PFAS recoveries were within 70-130%.
- The SW and WW matrices resulted in slightly low recoveries for long-chain PFAS.
 - PFUnA
 - 11CI-PF3OUdS
 - PFDoA
- PFMPA recoveries depended on the selected IDA.



EPA 533 Analyte %RSD




- Single lab RSD were $\leq 20\%$.
- The multi-lab PFAS RSD were generally $\leq 20\%$.
- 8:2 FTS = 22% for SW2
- PFBA spiked at 2 ng/L = 28% for RW1
 - Lab D MRL = 5 ng/L



EPA 537.1 Surrogate Recoveries and %RSD



- The multi-lab SS recoveries were within 70-130%.
- The SS recoveries slightly trended lower from PFHxS-13C2 to NETFOSAA-d5.
- The multi-lab RSD were within 10%.




EPA 537.1 Analyte Recoveries



- The multi-lab PFAS recoveries were within 70-130% except
 - 11Cl-PF3OUdS
 - PFDoA
 - PFTTrDA
 - PFTeDA


- SW1 and SW2 resulted in lower recoveries for these long-chain PFAS.



EPA 537.1 Analyte %RSD



- Single lab RSD were <20%.
- The multi-lab PFAS RSD were generally <20% except
 - PFDoA
 - PFTTrDA
 - PFTeDA
- SW2 resulted in higher RSD for these long-chain PFAS.
 - SW2 contained the highest TOC and HPC contents.




EPA 537.1 Experience



- **Low recoveries (50-69%) of SS NEtFOSAA-d5 were often observed for a small fraction of field samples.**
 - **Reanalysis**
 - **Re-extraction**
 - **Re-collection**


- **Long-chain PFAS recoveries could be improved by enhancing elution.**



Bottle Surface Adsorption Losses




8-oz HDPP bottles spiked at 400 ng/L
(Unpreserved, no chlorine, and stored at 1-6 °C)



Bottle Surface Adsorption Losses (Cont'd)



8-oz HDPE bottles spiked at 400 ng/L
(Unpreserved, no chlorine, stored at 1-6 °C)



Conclusions



- **EPA 537.1 and 533 are sensitive and robust drinking water methods.**
 - **Batch QC failures are rare.**
 - **2 ng/L or lower MRLs can be achieved for most PFAS of interest. Slightly higher MRLs (PFBA, PFHpA, 6:2 FTS, etc.) may be applicable for some labs.**
 - **Field reagent blank (FRB) contamination was rarely observed. Most FRB failures were related to switching bottles.**
 - **Carryover contamination from very high PFAS concentration samples was occasionally observed.**
 - **The most common EPA 537.1 QC failure was relatively low SS-NEtFOSAA-d5 recoveries, < 70%.**

Conclusions (Cont'd)



- EPA 537.1 and 533 are generally applicable for source water assessments – Pristine GW, SW, and treated WW matrices.
- EPA 533 is more robust in tolerating matrix interferences.
 - ESI suppression to PFBA and its IPS and IDA can be caused by high concentration common inorganic anions and/or polar organic compounds extracted particularly from WWTP effluents.
 - Sufficient rinsing of sample bottles and SPE cartridges with reagent water may be necessary.
 - It may be critical to choose appropriate IDA for those PFAS without their own labeled analogues available (e.g., PFMPA).

Conclusions (Cont'd)



- **EPA 537.1 is less robust in extracting long-chain PFAS (C10 and above) from difficult SW and WW matrices.**
 - **Adsorption losses of long-chain PFAS are a primary challenge, which can be enhanced by high concentration TOC and biological contents in SW and WW matrices.**
 - **Enhanced rinsing of sample bottles and SPE cartridge elution may be necessary.**
 - **Improved recoveries of long-chain PFAS was observed by increasing the elution solvent volumes.**
- **Most contents will be published in AWWA Water Science.**



Yongtao (Bruce) Li
Yongtao.Li@EurofinsET.com
574.472.5562
Eurofins Eaton Analytical, LLC
www.EurofinsUS.com