

NEMC, August 5, 2025

Chloronitramide Anion: A New DBP Detected in Drinking Water

Yongtao (Bruce) Li
Eurofins Eaton Analytical, LLC



Environment Testing

Acknowledgments



Joshua Whitaker

Dr. Julian L. Fairey

Jason Thornhill

Dr. David L. Wahman

Dr. Juliana R. Laszakovits

Eurofins Eaton Analytical

Department of Civil Engineering
University of Arkansas

Department of Civil Engineering
University of Arkansas

USEPA

Institute of Biogeochemistry and Pollutant
Dynamics, Switzerland

Presentation Outline



- Introduction
- Method Development and Validation
- Matrix Interferences
- Field Sample Analysis and Results
- Conclusions

Haloacetic acids (HAA5) MCL = 60 µg/L
Total trihalomethanes (TTHMs) MCL = 80 µg/L
Bromate MCL = 10 µg/L
Chlorite MCL = 1.0 mg/L
Chloramines as Cl₂ MRDL = 4.0 mg/L
Chlorine as Cl₂ MRDL = 4.0 mg/L
Chlorine dioxide as ClO₂ MRDL = 0.8 mg/L

Unregulated haloacetic acids & halomethanes
Haloacetonitriles & halonitromethanes
Haloketones & haloamides
Halofuranones such as MX
Aldehydes & ketones
Nitrosamines, etc.



Drinking Water D/DBPs

Unregulated Oxyhalides
Perchlorate
Chlorate
Iodate
Cyanogen chloride, etc.

Halocyclopentadienes
Chloronitramide anion
?

What is chloronitramide anion?

- Chloronitramide anion is a previously-unknown chemical compound in drinking water, which was recently identified as a decomposition product of inorganic chloramines by Fairey, et al. *Science*, 386(6724):882-887, (2024). <https://doi.org/10.1126/science.adk6749>


Name	Chloronitramide anion	Chloronitramide
Chemical formula	ClN_2O_2^-	ClHN_2O_2
Molecular weight (MW)	95.46 g/mol	96.47 g/mol
CAS No.	None	58999-86-3
Chemical structure		
Chemical property	Oxidant	Oxidant

Concern about chloronitramide anion?

- It was detected in all US municipal systems using chloramines to kill bacteria and viruses. The concentrations were greater than the MCLs of HAAs and TTHMs in a significant fraction of the samples studied.
- While the toxicity of chloronitramide anion is currently unknown, its structural similarity to other potentially toxic compounds raise concerns about potential health risks.
- The discovery of this previously-unknown DBP highlights the need for ongoing research into the safety of water treatment methods and the potential for unintended consequences.
- Researchers are calling for immediate toxicological evaluation, new analytical method development, and quantification of chloronitramide anion in source waters, finished drinking waters, and wastewater effluents as well.
- Future research may also advance in modeling water treatment and understanding DBP formation from the use of chloramines in drinking water disinfection.

Suitable Analytical Techniques




- Ion chromatography (IC) with electrical conductivity detection (ECD) or UV detection at 243 nm.
 - IC–electrospray ionization–ultrahigh resolution mass spectrometry (IC–ESI–UHRMS).
 - Hydrophilic interaction liquid chromatography–ultrahigh-resolution tandem mass spectrometry (HILIC–UHRMS).
 - Hydrophilic interaction liquid chromatography–tandem mass spectrometry (HILIC–MS/MS).
 - Reversed-phase liquid chromatography–tandem mass spectrometry (RPLC–MS/MS).
- 

Method Development & Validation

- RPLC–MS/MS
- Direct aqueous injection (DAI)
- Internal standard calibration
- Dark color containers & vials
- No sample preservation
- Fast analysis

DAI-RPLC-MS/MS Process





Chromatograms

m/z 95 > 51
 $^{35}\text{Cl}^{14}\text{N}_2^{16}\text{O}_2^- > ^{35}\text{Cl}^{16}\text{O}$

m/z 97 > 36.9
 $^{37}\text{Cl}^{14}\text{N}_2^{16}\text{O}_2^- > ^{37}\text{Cl}$

m/z 95 > 34.9
 $^{35}\text{Cl}^{14}\text{N}_2^{16}\text{O}_2^- > ^{35}\text{Cl}$

m/z 97 > 34.9
 $^{35}\text{Cl}^{15}\text{N}_2^{16}\text{O}_2^- > ^{35}\text{Cl}$



$\text{HCO}_3^- > \text{HCO}_3^-$

m/z 61 > 61

$\text{HSO}_4^- > \text{HSO}_4^-$

m/z 97 > 97

$\text{Cl}^- > \text{Cl}^-$

m/z 35 > 35

$\text{H}_2\text{PO}_4^- > \text{H}_2\text{PO}_4^-$

m/z 97 > 97

$\text{NO}_3^- > \text{NO}_3^-$

m/z 62 > 62

IDC – Sensitivity

Analyte Mass (m/z)	Spike Conc (µg/L) n ≥ 4	LCMRL (µg/L)	DL (µg/L)	MDL (µg/L)
95>34.9	0, 0.025, 0.005, 0.01, 0.02, 0.03, 0.04, 0.05, 0.1	0.0085	0.0030	0.0069
97>36.9	0, 0.005, 0.01, 0.02, 0.03, 0.04, 0.05, 0.1	0.0079	0.0045	0.0096
95>51	0, 0.04, 0.05, 0.1, 0.2, 0.3, 0.4, 0.5	0.049	0.0076	0.049

LCMRL = Lowest Concentration Minimum Reporting Level

MRL Confirmation (n = 7)

Analyte Mass (m/z)	Spike Conc (µg/L)	Upper PIR (%)	Lower PIR (%)	Mean Rec (%)	RSD (%)
95>34.9	0.005	130 ≤ 150	83 ≥ 50	106.2	5.6
97>36.9	0.01	133 ≤ 150	58 ≥ 50	96.0	9.9
95>51	0.05	145 ≤ 150	54 ≥ 50	99.4	11.6

PIR = Prediction Interval of Results

IDC – Accuracy & Precision

QC Sample	Spike Conc (µg/L)	Analyte/IS Mean Rec ± RSD (%), n = 7			
		Mass m/z 95>34.9	Mass m/z 97>36.9	Mass m/z 95>51	* IS Mass m/z 97>34.9
LCS	25.0	100.6 ± 0.6	100.3 ± 0.3	100.7 ± 0.9	98.3 ± 0.6
LCS	5.0	100.2 ± 1.1	100.9 ± 1.1	104.5 ± 2.0	92.0 ± 4.3
LCS	0.1	100.5 ± 1.7	104.4 ± 0.8	92.8 ± 7.0	97.2 ± 0.4

Impacts of pH Values

DI-Water pH Value	Spike Conc (µg/L)	Analyte/IS Mean Rec ± RSD (%), n = 4			
		Mass m/z 95>34.9	Mass m/z 97>36.9	Mass m/z 95>51	* IS Mass m/z 97>34.9
pH 10	5.0	100.1 ± 3.0	100.8 ± 3.1	100.8 ± 3.0	78.9 ± 3.0
pH 9	5.0	98.9 ± 4.8	99.7 ± 4.7	100.3 ± 4.6	79.7 ± 4.7
pH 8	5.0	108.3 ± 11.2	109.0 ± 11.0	109.3 ± 11.5	79.3 ± 2.4
pH 7	5.0	102.0 ± 4.1	102.6 ± 4.2	103.4 ± 4.4	79.6 ± 1.3
pH 6	5.0	102.1 ± 1.8	102.9 ± 1.8	103.3 ± 2.0	78.2 ± 2.3
pH 5	5.0	106.8 ± 3.8	107.3 ± 3.8	107.4 ± 3.6	75.0 ± 3.2
pH 4	5.0	103.5 ± 1.5	103.9 ± 1.8	104.6 ± 1.7	80.1 ± 1.5

Impacts of Free Chlorine Levels

Free Chlorine Level (as Cl ₂)	Spike Conc (µg/L)	Analyte/IS Mean Rec ± RSD (%), n = 4			
		Mass m/z 95>34.9	Mass m/z 97>36.9	Mass m/z 95>51	* IS Mass m/z 97>34.9
0.5 mg/L	5.0	97.6 ± 0.7	97.2 ± 0.6	98.6 ± 0.8	108.2 ± 0.9
1.0 mg/L	5.0	97.9 ± 0.2	97.5 ± 0.2	98.4 ± 0.4	105.6 ± 0.3
2.0 mg/L	5.0	97.1 ± 0.7	97.1 ± 0.6	98.1 ± 0.9	103.6 ± 0.1
3.0 mg/L	5.0	96.3 ± 0.5	95.9 ± 0.6	97.3 ± 0.7	101.9 ± 0.3
4.0 mg/L	5.0	96.4 ± 0.9	95.8 ± 0.6	97.0 ± 1.0	99.8 ± 0.4

Impacts of Chloramine Concentrations

4:1 Cl ₂ :N Concentration	Analyte/IS Mean Conc ± STD (µg/L), n = 4			
	Mass m/z 95>34.9	Mass m/z 97>36.9	Mass m/z 95>51	* IS Mass m/z 97>34.9
0.5 ppm Cl ₂ :0.125 ppm N	1.20 ± 0.01	1.22 ± 0.02	1.22 ± 0.01	10.0 ± 0.2
1 ppm Cl ₂ :0.25 ppm N	2.40 ± 0.06	2.45 ± 0.06	2.48 ± 0.07	9.8 ± 0.2
2 ppm Cl ₂ :0.5 ppm N	2.15 ± 0.78	2.15 ± 0.79	2.21 ± 0.80	10.1 ± 0.2
4 ppm Cl ₂ :1 ppm N	3.88 ± 1.16	3.80 ± 1.20	3.94 ± 1.18	9.9 ± 0.1
6 ppm Cl ₂ :1.5 ppm N	6.88 ± 0.24	6.86 ± 0.25	6.96 ± 0.26	9.6 ± 0.3
8 ppm Cl ₂ :2 ppm N	9.22 ± 0.03	9.16 ± 0.03	9.28 ± 0.03	9.7 ± 0.1
10 ppm Cl ₂ :2.5 ppm N	9.37 ± 0.04	9.37 ± 0.08	9.46 ± 0.05	9.6 ± 0.1

Impacts of Common Anions

DI-Water pH Value	Spike Conc (µg/L)	Analyte/IS Mean Rec ± RSD (%), n = 4 or 7			
		Mass m/z 95>34.9	Mass m/z 97>36.9	Mass m/z 95>51	* IS Mass m/z 97>34.9
Chloride at 250 mg/L	5.0	101.4 ± 1.5	101.5 ± 0.7	103.4 ± 0.8	63.8 ± 1.3
Sulfate at 250 mg/L	5.0	93.3 ± 0.7	93.5 ± 0.6	95.5 ± 0.6	78.3 ± 0.3
Nitrate at 10 mg N/L	5.0	92.5 ± 0.4	93.2 ± 0.5	94.5 ± 0.5	72.6 ± 0.2
Phosphate at 50 mg P/L	5.0	95.8 ± 0.2	96.1 ± 0.3	97.3 ± 0.5	78.3 ± 0.3
Bicarbonate at 400 mg/L	5.0	95.8 ± 3.1	97.9 ± 2.8	98.4 ± 2.4	7.8 ± 2.4
LFSSM	25.0	94.8 ± 0.2	94.3 ± 0.3	95.9 ± 0.2	99.6 ± 0.6
LFSSM	5.0	95.9 ± 0.4	95.9 ± 0.4	104.9 ± 0.5	104.7 ± 0.5
LFSSM	0.1	99.9 ± 6.1	105.2 ± 4.5	98.6 ± 19.4	103.8 ± 1.3

Initial Results of Cl₂/UV/H₂O₂ DW (1–6 °C)

Sample Site	Holding Time (Day)	Chloronitramide Anion Mean Conc of FS & DUP (µg/L)		
		Mass m/z 95>34.9	Mass m/z 97>36.9	Mass m/z 95>51
AZ-DS-1	3	2.26	2.28	2.33
AZ-DS-2	7	0.56	0.52	0.56
CA-DS-10	11	0.32	0.28	0.33
CA-DS-11	11	0.16	0.16	0.14
CA-DS-9	2	0.24	0.25	0.26
CO-DS-2	3	0.41	0.41	0.41
IN-DS-5	6	32.6	32.7	32.9
IN-DS-6	0	0.85	0.86	0.86
IN-DS-7	2	1.07	1.08	1.07
IN-DS-8	0	0.81	0.82	0.81
IN-DS-9	5	2.15	2.16	2.23
MI-DS-1	5	<0.1	<0.1	<0.1

Initial Results of Chloraminated DW (1–6 °C)

Sample Site	Holding Time (Day)	Chloronitramide Anion Mean Conc of FS & DUP (µg/L)		
		Mass m/z 95>34.9	Mass m/z 97>36.9	Mass m/z 95>51
CA-DS-1	1	111	111	112
CA-DS-2	1	28.2	28.3	28.5
CA-DS-3	7	19.5	19.52	19.8
CA-DS-4	7	22.1	22.2	22.5
CA-DS-5	6	28.2	28.3	28.7
CA-DS-6	6	21.8	21.8	22.1
CA-DS-7	6	11.8	11.9	12.2
CA-DS-8	15	0.85	0.85	0.83
CO-DS-1	2	28.1	28.2	28.3
FL-DS-1	1	55.4	55.4	55.6
GA-DS-1	4	1.10	1.08	1.17
IN-DS-1	6	31.7	31.8	32.1
IN-DS-2	4	39.7	39.7	40.1


Initial Results of Chloraminated DW (1–6 °C)

Sample Site	Holding Time (Day)	Chloronitramide Anion Mean Conc of FS & DUP (µg/L)		
		Mass m/z 95>34.9	Mass m/z 97>36.9	Mass m/z 95>51
IN-DS-3	6	3.41	3.54	3.47
IN-DS-4	3	7.65	7.72	7.77
KY-EP-1	2	5.34	5.35	5.46
KY-DS-1	2	8.12	8.17	8.29
KY-EP-2	2	3.59	3.59	3.61
KY-DS-2	2	7.02	7.03	7.17
MA-DS-1	4	18.4	18.5	18.5
MN-DS-1	2	12.0	12.0	12.3
MN-DS-2	2	7.67	7.75	7.94
MO-DS-1	7	14.2	14.1	14.7
TX-DS-1	3	91.0	90.5	92.0
TX-DS-2	7	83.1	82.5	82.9
TX-DS-3	5	145	145	146


MS Recoveries (stored at 1-6 °C)

Sample Site	Native Conc (µg/L)	Spike Conc (µg/L)	Analyte Mean Rec ± RSD (%), n = 4		
			Mass m/z 95>34.9	Mass m/z 97>36.9	Mass m/z 95>51
AZ-DS-1	2.13	5.0	103.3 ± 1.2	104.0 ± 1.0	104.7 ± 1.8
CA-DS-1	122	100	98.5 ± 1.2	97.0 ± 1.0	99.8 ± 1.0
CA-DS-4	41.6	20.0	105.5 ± 0.8	104.7 ± 0.9	105.6 ± 0.7
CO-DS-2	0.51	5.0	97.4 ± 0.8	97.8 ± 0.6	98.0 ± 1.0
IN-DS-2	68.2	50.0	104.5 ± 2.7	103.5 ± 2.8	104.2 ± 2.7
MN-DS-1	21.2	10.0	99.6 ± 0.6	98.8 ± 0.6	97.6 ± 0.6


Impacts of Holding Time (1-6 °C)





Impacts of Storage Temperature (1)



Impacts of Storage Temperature (2)



Impacts of Storage Temperature (3)



MS Recoveries (stored at -25 °C)

Sample Site	Native Conc (µg/L)	Spike Conc (µg/L)	Mean MS/MSD Conc (ug/L)	Mean MS/MSD Rec (%)	MS/MSD RPD (%)
AZ-DS-1	2.26	5	8.1	117.3	0.7
CO-DS-2	0.41	5	5.6	104.4	1.5
FL-DS-1	52.8	25	76.6	95.2	4.5
CA-DS-1	111	100	224	112.6	1.8
CA-DS-4	77.3	100	168	90.6	2.1
CA-DS-6	46.4	100	139	92.8	0.2
IN-DS-1	89.5	100	192	102.0	1.3
IN-DS-2	88.3	100	168	79.4	2.4
MN-DS-1	200	100	297	96.4	4.7

Conclusions

- The DAI–LC–MS/MS method could determine chloronitramide anion in a typical alkalinity range of drinking water at an MRL of 0.1 µg/L.
- This method demonstrated satisfactory sensitivity, accuracy and precision. Common anions except bicarbonate, pH, free chlorine concentrations, and chloramine concentrations did not interfere with the quantitation of chloronitramide anion.
- Surprisingly, freezing samples caused chloramine decomposition and more chloronitramide anion formation.
- Chlorine/chloramine quenching agents could not be used to preserve drinking water samples. Storing samples in a freezer was not recommended.

Conclusions (Cont'd)

- Chloronitramide anion was detected in all studied chloraminated drinking water samples at the initial analysis concentrations of 0.85 to 111 µg/L.
- The concentrations substantially increased (up to 200 µg/L) with the increase in holding time due to the continuous decomposition of chloramines. Ideally, the samples should be analyzed as early as possible, within 7 days of holding time.
- 11 out of 12 chlorinated drinking water samples were also detected with chloronitramide anion at low µg/L with an exception for MI-DS-1.
- Chloronitramide anion was initially detected in chlorinated IN-DS-5 at 32.6 µg/L, which was increased to ~90 µg/L after two weeks of holding time or freeze/thaw.

Yongtao (Bruce) Li

Yongtao.Li@ET.EurofinsUS.com

574.472.5562

THANK YOU



Environment Testing