

Ion Chromatography Based Isolation and Quantitation of Chloronitramide Anion in Drinking Water

Jason A. Thornhill and Julian L. Fairey

Department of Civil Engineering, University of Arkansas, Fayetteville, AR 72701

National Environmental Monitoring Conference, August 4-7, 2025, St. Louis, MO


Chloronitramide anion (Cl-N-NO_2^-) is a recently identified decomposition product of inorganic chloramines

- Chloramines are used to disinfect tap water of ~113 M in the US alone
- Cl-N-NO_2^- forms at concentrations up to ~130 $\mu\text{g}\cdot\text{L}^{-1}$
- Current analytical method for quantitation uses HILIC–UHRMS (Hydrophilic Interaction Liquid Chromatography–Ultrahigh Resolution Tandem Mass Spectrometry)
- Need simpler quantitation method for broad occurrence studies




Incomplete Cl–N–NO₂[–] Isolation at 0.5 mL•min^{–1}

- Two flowrates
 - 0.5 and 0.7 mL•min^{–1}
- A: Baseline does not reach zero
 - Incomplete Cl–N–NO₂[–] separation from common anions
- B: UV₂₅₄ shows increasing elution time with decreasing flowrate by ~30 min
 - 0.7 mL•min^{–1} selected for isolation




Incomplete Cl–N–NO₂[–] Isolation with 1 mL Injection Loop

- Three eluent strengths
 - 1000X, 500X, 2000X
- B: Baseline does not reach zero
 - Incomplete Cl–N–NO₂[–] separation from common anions
- C & D: UV₂₅₄ shows increasing elution time with decreasing eluent strength by ~60 min
 - 1000X selected for isolation




Complete Cl-N-NO₂⁻ Isolation with 100 μL Injection Loop

- C: Cl-N-NO₂⁻ elutes ~40 min with 100 μL injection loop
- B: Baseline reaches zero prior to 40 min with 100 μL injection loop
 - Complete separation
- B: 1 mL injection loop overloading column
 - 100 μL injection loop selected for isolation




Decreased Elution Time with Increased Column Temperature

- B: Cl-N-NO_2^- elutes ~32 & 38 min, depending on concentration
- B: **Peak** height increases with increasing column temperature
- B: Elution time decreases with increasing column temperature
- Higher temperature increases IC equilibration time by 30–60 min
- 50 °C selected for isolation



Collection Interval of Cl-N-NO_2^-

- A: IC–EC
 - Complete isolation from common anions
 - Start = 34 min
 - End = 45 min
- B: IC–UV₂₄₃
 - Collection start at peak front
 - Collection end after signal < 50 mAU




Cl–N–NO₂[–] Isolation: IC Method

Isolation		
Model	Metrohm 850 Professional IC	
Injection Loop	100	μL
Column	Metrosep A Supp 7 – 250/4.0	
Temperature	50	°C
Flow Rate	0.7	mL•min ^{–1}
Eluent	Metrohm A Supp 7 Eluent	
Eluent Strength	1000X	
Detectors	Metrohm IC Conductivity Detector	
	Metrohm 887 Professional UV/Vis Detector	
UV Wavelength	243	nm
UV Bandwidth	2	nm


Cl-N-NO₂⁻ Quantitation: Minimal Impact of Column Temperature

- A & C: Cl-N-NO₂⁻ elution ~50 min with 1 mL injection loop
 - 1 mL injection loop used to maximize detection
- B & D: Cl-N-NO₂⁻ eluted ~50 min at 45 and 55 °C
- **Column** temperature does not impact elution time or peak height




Cl-N-NO₂⁻ Elution Decreases with Increasing Eluent Strength

- A & C: Cl-N-NO₂⁻ elutes ~75 min with 1000X eluent and ~55 min with 500X eluent
- B & D: 500X eluent selected to increase throughput during quantitation




UV₂₄₃ Selected for Cl–N–NO₂[–] Quantitation

- A & C: Cl–N–NO₂[–] elutes ~45 min
- B: No variation in peak height with triplicate injection
- C: UV₂₀₀ excluded due to detection of other constituents (nitrate)
- D: Cl–N–NO₂[–] UV molar absorptivity maxima at 243 nm selected for quantitation because S/N > 254 nm



Cl-N-NO₂⁻ Peak Tailing Necessitated Peak Area for Quantitation

- A & C: Cl-N-NO₂⁻ elutes ~43 min
- B & D: As concentration increases, peak tailing becomes more apparent
- B: IC-EC area taken 0.5 min prior to and following due to flat, noise-free baseline
- D: IC-UV₂₄₃ taken at start and end of peak due to baseline noise



Cl-N-NO₂⁻ Quantitation: IC Method

Quantitation		
Model	Metrohm 850 Professional IC	
Injection Loop	1	mL
Column	Metrosep A Supp 7 – 250/4.0	
Temperature	45	°C
Flow Rate	0.7	mL•min ⁻¹
Eluent	Metrohm A Supp 7 Eluent	
Eluent Strength	500X	
Detectors	Metrohm IC Conductivity Detector	
	Metrohm 887 Professional UV/Vis Detector	
UV Wavelength	243	nm
UV Bandwidth	2	nm


Cl-N-NO₂⁻ Standard Curve by IC-EC and IC-UV₂₄₃

A. IC-EC

- Standards range from 20–300 µg•L⁻¹

B. IC-UV₂₄₃

- Standards range from 10–300 µg•L⁻¹
- Ten standard curves were made and those with R² > 0.990 averaged and used to determine concentrations




Cl-N-NO₂⁻ MDL, LOD, and LOQ

- MDL determined from EPA method
- IC-UV₂₄₃ LOD and LOQ from blank determination
- IC-EC LOD and LOQ from linear regression
- IC based LOD and LOQ > HILIC-UHRMS
- IC-UV₂₄₃ has lower MDL and LOQ than IC-EC

	MDL (μg•L ⁻¹)	LOD (μg•L ⁻¹)	LOQ (μg•L ⁻¹)
IC-EC	22.4	11.8	39.4
IC-UV ₂₄₃	13.0	10.4	12.7
HILIC-UHRMS		0.17	0.58


IC-UV₂₄₃ has Lower Reproducibility Variation for Cl-N-NO₂⁻

- Testing over 6 days
- A & B: IC-EC quantitation accurate at 25 & 50 $\mu\text{g}\cdot\text{L}^{-1}$ but had negative residuals at 100 & 200 $\mu\text{g}\cdot\text{L}^{-1}$ (< ~20 %)
- C & D: IC-UV₂₄₃ quantitation like IC-EC, although residuals were lower at high concentrations (< ~15 %)



Cl-N-NO₂⁻ Quantitation in Tap Waters and Synthetic Waters

- A: Comparison between IC methods and HILIC-UHRMS
- B: Residuals calculated from IC minus HILIC-UHRMS
- Black line is 1:1 line and dashed lines are $\pm 20 \mu\text{g}\cdot\text{L}^{-1}$
- IC-EC residuals are closer to zero compared to IC-UV₂₄₃



Conclusions and Future Work

Conclusions:

- IC–EC LOQ = $39.4 \mu\text{g}\cdot\text{L}^{-1}$ with reproducibility of 20% for lab-grade waters and within $\sim 10\%$ of HILIC–UHRMS
- IC–UV₂₄₃ LOQ = $12.7 \mu\text{g}\cdot\text{L}^{-1}$ with reproducibility of 15% for lab-grade waters and within $\sim 50\%$ of HILIC–UHRMS

Future Work:

- Repeatability testing with water matrices spiked with Cl–N–NO₂[–] to determine matrix effects
- Formation conditions for Cl–N–NO₂[–] in drinking water systems